

Autoclave Gloves

Extremely tough & hard wearing whilst offering comfort and flexibility.



Product Code: ACG 001

Our Autoclave Kevlar gloves are tough, durable, comfortable, flexible, water resistant and capable of withstanding temperatures up to 450°C.

The glove's facing material is Kevlar backed with aluminised para-aramid fibre at 500g/m².

The glove seams are stitched with Kevlar wire and have an internal fleece insulation in wool or cotton.

The gloves contain no harmful materials or substances and comply with EC type-examination certificate No 0072/217/162/07/95/0378 / ext 01 and 02 Standard EN 420.

Available in sizes 7,8,10,11

