

## Distal Soak

**Clear visual instrument tip protection.**



**Product Code: 50058, 50065, 50069**

The Distal Soak is a highly effective method of ensuring that the delicate tips of surgical instruments remain protected and moist after use.

By preventing the drying out of the tip, any contamination becomes subsequently easier to remove, making the cleaning and sterilisation processes more effective.

The vials are available as closed lid - allowing the vials to be pre-charged with sterile solution, or open lid - allowing vials to be filled once the distal tip has been inserted.

The vials are fashioned from polypropylene, with silicone lids and a Plexiglas® tube holder (sold separately) is available for proper handling.

