

Face Shield

Light weight facial protection providing comprehensive cover and ventilation



Product Code: 18000W-100, 18000-100, 16000-100

The Minerva range of face shields are ideal for use in departments where maximum personal protection is the goal, providing the best splash protection against exposure to blood borne pathogens.

The standard face shield measures 9 inches high and 13 inches wide and is also available with a wider fit, with or without an additional drape.

The full face shield with drape adds a 13 inch fluid barrier that extends below the integral full face shield. The barrier fabric tucks inside a sterile gown, providing additional protection to prevent blood from splashing under the shield.

Air circulates around the fabric to allow for easy breathing and preventing moisture accumulation that can lead to fogging.

The shields are CE marked and manufactured under ISO conditions for the highest quality & delivers OSHA level protection against face splash.