


Instrument Tip Protectors

Pocketed, channelled, vented and flexible tip covers for maximum protection



The Famos Instrument Protector provides a secure cover for sharps.

They measure 50 x 125 mm and are designed for sterilisation conditions.



Vented instrument guards are designed to accommodate full impact sterilisation yet provide maximum tip protection.

The see-through covers means tips are easily visible with a design that prevents the protector from sticking to the instrument.

The channel design tip protectors prevent moisture from being trapped between the protector and the instrument.

Being vented, they ensure that moisture quickly evaporates during the sterilisation cycle.



Mesh style tip protectors (sleeved or closed tip) are designed for instruments of a more delicate nature or design.

These lightweight durable polyethylene single use tip protectors are perfect for instruments such as rigid endoscopes.