

Stainless Steel Spray Gun

A sleek design matched by tough, hard wearing reliable performance



Product Code: 24.901.080

This spray gun is constructed entirely of stainless steel making it one of the most durable spray guns on the market.

It's design and construction makes it suitable for use with demineralised water as well as air and being capable of being entirely sterilised ready for use in the operating theatre or the clean room.

The stainless steel spray gun comes equipped with two nozzles, an injection attachment piece, a trigger locking ring, a wall-mount, a hose with stainless steel snap-closing connector and a coupling protection cover.

Additional optional attachments are also available