

The Standard Set Spray Gun

A spray gun for water or air with smooth flow control across the entire trigger range



Product Code: 24.900.350 (Spray gun)

The standard set spray gun by Famos is the ideal instrument for cleaning laboratory glassware, pipettes, syringes, catheters, cuvettes, endoscopes, cystoscopes, cannulae and a variety of other instruments and articles used in clinics, laboratories and medical practices.

The spray gun features a new type of internal flow control mechanism which permits the smooth and simple adjustment of the water or air flow across the entire range of the trigger movement.

The desired maximum water or airflow can be set by using the small knurled screw located behind the trigger.

A selection of optional attachments are available & are connected by seating them on the safety cone located at the tip.

The attachment fit securely without drama & high pressure cannot disconnect them.

Equipped with a heat insulated grip, the spray gun is incredibly easy to use.

